


TD-9030 Bluetooth Dongle Operation Instruction

Ver 1.0 2015-12 311-9030300-001

Dongle Overview

- 1 Cap
- 2 Strip Port Connector
- 3 Bluetooth Indicator



Dongle Operation

1. Start with the meter off. Take off the cap and insert the Bluetooth dongle into a strip port on the meter. Make sure the Bluetooth indicator is on the right side up before insertion.
2. The Bluetooth indicator on the dongle will be in flashing blue and "PC" will appear on the meter display, indicating that the Bluetooth function is on.
3. And then the Bluetooth light will turn to solid blue, indicating that the Bluetooth connection is established.

Bluetooth Indicator on the Dongle

BLUETOOTH INDICATOR	STATUS
Flashing blue	The Bluetooth function is on and waiting for connection.
Solid blue	The Bluetooth connection is established.

Symbol Information

SYMBOL	REFERENT	SYMBOL	REFERENT
	Consult instructions for use		Humidity Limitation
	Temperature limitation		Collection for electrical and electronic equipment
SN	Serial number	CE	CE mark
	Caution, consult accompanying documents		Manufacturer

Specification

Model No.	TD-9030
Dimension & Weight	43.2 (L) x 15 (W) x 5 (H) mm, 3.56 g
Data Interface	Strip port
External Output	Bluetooth
Power Source	NA (provided from meter)
Frequency	2.4 GHz
Operating Condition	0°C to 50°C (32°F to 122°F), 15% to 93% R.H. (non-condensing)
Storage Condition	-20°C to 60°C (-4°F to 140°F), 10% to 95% R.H. (non-condensing)

TaiDoc Technology Corporation
B1-7F, No. 127, Wugong 2nd Rd., Wugu Dist.,
24888 New Taipei City, Taiwan
www.taidoc.com

